

Accessory compatibility guide

Effective date: September 2018

| Product | Tempus platform | North platform | Era platform |
|---|-----------------|----------------|--------------|
| For patients: | | | |
| Rechargeable kit¹ | • | • | |
| Remote Plus app² | • | | |
| uControl app | • | • | |
| TV Connector² | • | | |
| Remote control 2 | • | • | |
| uStream | • | • | |
| uDirect 3 | • | • | |
| uTV 3³ | • | • | • |
| uMic³ | • | • | • |
| IIC remote (micro CIC remote) | • | • | • |
| Remote control | | | • |
| Smart Control remote | | | • |
| uDirect 2 | | | • |
| uTV 2³ | | | • |
| For hearing healthcare professionals: | | | |
| iCube II programming | • | • | • |
| iCube programming | | | • |

¹Product-specific rechargeable kit offered for Max SP R, Moxi All R, Moxi Fit R and Stride M R hearing instruments

²Compatible with Moxi All, 2.4 GHz hearing instruments

³uStream, uDirect 3 or uDirect 2 required

uControl app, uStream, uDirect 3, uTV 3, uMic, and IIC remote (micro CIC remote) are not compatible with the Moxi All hearing instruments
Accessory compatibility may also be limited by product style and technology level.

| Tempus platform: | | North platform: | | Era platform: | |
|------------------|-----------------|-----------------|--------------------------|--------------------------------|--|
| T Max SP | T Stride P Dura | Stride M | Moxi Dura | Moxi ² | Quantum ² wireless ITEs |
| T Max UP | T Moxi Now | Stride P | Stride wireless ITEs | Moxi ² Kiss | Quantum ² non-wireless ITEs |
| T Moxi All | T Moxi Kiss | Stride P Dura | Stride non-wireless ITEs | Moxi ² Dura | Max SPm |
| T Insera | T Moxi Fit | Moxi Now | Stride IIC | Quantum ² S | Max SP |
| T Stride M | T Moxi Dura | Moxi Kiss | | Quantum ² HP | Shine Rev BTEs |
| T Stride P | | Moxi Fit | | Quantum ² micro CIC | Shine Rev ITEs |

Accessory pairing guide

| Product | uStream | uDirect 3 | uDirect 2 |
|---------------------------------|----------------------|-------------------------------|----------------------|
| Remote Plus app | | | |
| uControl app | Remote control | Remote control and My Ratings | No |
| uTV 3 | Streaming & charging | Streaming & charging | Streaming & charging |
| uMic | Streaming | Streaming | Streaming |
| uTV 2 | Streaming | Streaming & charging | Streaming & charging |