

Accessory compatibility guide

Product	Tempus platform	North platform	Era platform
Rechargeable kit ¹	•	•	
uControl app	•	•	
Remote control 2	•	•	
uStream	•	•	
uDirect 3	•	•	
uTV 3	• ²	• ²	• ²
uMic	• ²	• ²	• ²
IIC remote (micro CIC remote)	•	•	•
Remote control			•
Smart Control remote			•
uDirect 2			•
uTV 2			• ²
iCube II programming	•	•	•
iCube programming			•

¹compatible with Moxi Fit R and Stride M R hearing instruments on Tempus and North platforms

²uStream, uDirect 3 or uDirect 2 required.

Accessory compatibility may also be limited by product style and technology level.

Tempus platform:	North platform:		Era platform:	
T Insera	Stride M	Moxi Fit	Moxi ²	Quantum ² micro CIC
T Stride M	Stride P	Moxi Dura	Moxi ² Kiss	Max SPm
T Stride P	Stride P Dura	Stride wireless ITEs	Moxi ² Dura	Max SP
T Stride P Dura	Moxi Now	Stride non-wireless ITEs	Quantum ² S	Shine Rev BTEs
T Moxi Now	Moxi Kiss	Stride IIC	Quantum ² HP	Shine Rev ITEs
T Moxi Kiss			Quantum ² wireless ITEs	
T Moxi Fit			Quantum ² non-wireless ITEs	
T Moxi Dura				

Accessory pairing guide

Product	uStream	uDirect 3	uDirect 2
uControl app	Remote control	Remote control and My Ratings	No
uTV 3	Streaming & charging	Streaming & charging	Streaming & charging
uMic	Streaming	Streaming	Streaming
uTV 2	Streaming	Streaming & charging	Streaming & charging