





Earmold Option Reference Chart *Stride and Max BTEs*



Product	Skeleton	Semi-Skeleton	Full Shell	Half Shell	Canal Lock	Canal	CROS	Hollow Sleeve	Solid Sleeve
Stride™ M and Stride M R BTE	✓	✓	✓	✓	✓	✓	✓	✓	✓
Stride P BTE	✓	✓	✓	✓	✓	✓	✓	✓	✓
Stride P Dura BTE	✓	✓	✓	✓	✓	✓	✓	✓	✓
Max™ SPm BTE	✓	✓	✓	✓	✓	✓			
Max SP BTE	✓	✓	✓	✓	✓	✓			
Material									
Acrylic	✓	✓	✓	✓	✓	✓	✓	✓	✓
Silicone	✓	✓	✓	✓	✓	✓	✓		✓
Soft Silicone			✓	✓		✓			✓
Finish									
Gloss	✓	✓	✓	✓	✓	✓	✓	✓	✓
Satin	✓	✓	✓	✓	✓	✓	✓	✓	✓
Color									
Glitter/Solid/Swirl	Earmolds available in 15+ color options. Refer to Price book (028-5164-03) or Unitron's Innovative Product Family brochure (028-6248-03) for more information.								
Tubing									
Slim Tube	✓	✓	✓	✓	✓	✓	✓	✓	✓
Power Slim Tube (Max only)	✓	✓	✓	✓	✓	✓			
13 Regular	✓	✓	✓	✓	✓	✓	✓		
13 Heavy Wall	✓	✓	✓	✓	✓	✓	✓		
QuickSnap™ Regular	✓	✓	✓	✓	✓	✓	✓		
QuickSnap Heavy Wall	✓	✓	✓	✓	✓	✓	✓		
Venting									
No Vent	✓	✓	✓	✓	✓	✓		✓	✓
IntelliVent	✓	✓	✓	✓	✓	✓		✓	✓
Pressure Vent	✓	✓	✓	✓	✓	✓		✓	✓
Select-A-Vent (2.0 mm) (2.5mm) (3.0 mm)	✓	✓	✓	✓	✓	✓		✓	✓
IROS Vent	✓	✓	✓	✓	✓	✓		✓	✓
As Large As Possible	✓	✓	✓	✓	✓	✓		✓	✓



	 Hollow Sleeve	 Solid Sleeve	 cShell (hard)	 cShell (soft)
Product				
Moxi™ Now	✓	✓	✓	✓
Moxi Fit and Moxi Fit R	✓	✓	✓	✓
Moxi All and Moxi All R	✓	✓	✓	✓
Moxi Kiss	✓	✓	✓	✓
Moxi Dura	✓	✓	✓	✓
Receiver				
Standard (xS)	✓	✓	✓	✓
Power (xP)			✓	✓
Super Power plus (xSP Plus)			✓	
Material				
Acrylic	✓	✓	✓	
Silicone		✓		✓
Soft Silicone		✓		
Finish				
Gloss	✓	✓	✓	✓
Satin	✓	✓	✓	✓
Color				
Pink	✓	✓	✓	
Brown	✓	✓	✓	
Tan	✓	✓	✓	
Clear	✓	✓	✓	✓
Cocoa	✓	✓	✓	
Blue/Red	✓	✓	✓	
Translucent pink	✓	✓		
Translucent brown	✓	✓		
Venting				
No Vent	✓	✓	✓	✓
IntelliVent	✓	✓	✓	✓
Pressure Vent	✓	✓	✓	✓
Select-A-Vent (2.0 mm) (2.5mm) (3.0 mm)	✓	✓	✓	✓
IROS Vent	✓	✓	✓	✓
As Large As Possible	✓	✓	✓	✓